

**X-Area**  
SINGLE CRYSTAL  
DIFFRACTION  
SOFTWARE


## THE STOE SINGLE CRYSTAL DIFFRACTION SOFTWARE PACKAGE



### SINGLE CRYSTAL DIFFRACTOMETRY

X-Area is the complete, easy to use, yet flexible and powerful area detector software package for data collection and data evaluation on the STOE STADIVARI and STOE IPDS diffractometer series

- Highest quality of unit cell and intensity data
- Providing SHELX-compatible intensity data sets
- Full integration of STOE's FaceIt<sup>Video</sup> and X-Shape
- Liberal license policy: free software updates for 3 years, license for unlimited use in your department



YOUR PARTNER IN X-RAY DIFFRACTION

STOE & Cie GmbH | [WWW.STOE.COM](http://WWW.STOE.COM)

# SINGLE CRYSTAL DIFFRACTION



## DATA ACQUISITION


### INSTRUMENT CONTROL / MEASUREMENT / RUN OPTIMIZER

- Intuitive interface to the diffractometer, offering direct access to all functions
- Easy centering of the sample with the aid of the integrated Faceit<sup>Video</sup> function
- Automated data collection is started with just a few steps in the measuring program
- Frames with 32 bits per pixel facilitate high data accuracy
- User-friendly GUI optimizes the run strategy (e.g. in terms of high data completeness according to the crystal system and orientation), allowing faster and more dedicated collection of the data

## PRESENTATION


### FRAME GRAPHICS

- X-Area offers a comprehensive graphics program for the inspection of collected frames
- Interactive display options provide an easy-to-use way for checking the quality of the crystal to be investigated, e.g. unexpected splitting of reflections can be examined easily
- Diffraction patterns can be checked for anomalous effects, as diffuse scattering




### CRYSTAL MONITORING

Video-camera signal to enable user-friendly crystal centering or even face indexing, if needed



### RUN OPTIMIZER

Powerful tool for enabling time-optimized data collection



### REFLECTION INSPECTION

Sophisticated graphics software for a closer look on the frames


# IMAGE PROCESSING

## INDEX / CELL / REFINE

- Peak-finding routine scans frames rapidly in multi-processing mode
- Resulting peaks can be indexed, either automatically or based on a robust graphics method
- Unit-cell parameters are refined according to the crystal system
- Unlimited number of peaks for the refinement


## INTEGRATION

- Use of elliptical reflection profiles and  $\alpha_1/\alpha_2$ -splitting option for integration process
- Automated optimization of integration parameters
- Self-acting detection of reflection overlap
- Graphical control over the integration process
- Masking of user-specified detector areas possible
- Support of high-pressure cells by automatically calculated shading masks
- SHELX-compatible intensity data file by default
- XD-compatible file as an option




### ROBUST INDEXING

Graphics supported indexing enables full control of the crystal quality



### CONTROL BY GRAPHICS

During indexing, special effects as satellites may be detected easily




### REFLECTION INTEGRATION

Reliable integration procedure provides accurate intensity data set

# DATA ANALYSIS


- E-statistics plot in order to decide whether space group is centrosymmetric
- Analyzer for convenient determination of the correct Laue group
- Automated space-group determination
- Display of peaks in reciprocal-space viewer
- Difficult patterns can be inspected with the aid of layer representations in a reciprocal-space coordinate system, built from pixels of the collected frames
- Transformation of frame pixels into "powder diagrams"



### LAUE ANALYZER

Quick and easy check of the Laue group

# SINGLE CRYSTAL DIFFRACTION



## X-SHAPE

Automated crystal-shape optimization for numerical absorption correction


## CORRECTIONS

- Lp and air-absorption correction
- Correction of absorption of x-rays by crystal (numerical correction or intensity scaling based on spherical harmonics in conjunction with symmetry-related reflections)
- Automated version of STOE's X-Shape
- Inter-frame scaling, based on polynomials
- Correction of crystal decomposition
- Rejection of outliers

## EXTENDED ISSUES


### MULTI-DOMAIN SYSTEMS

- Semi-automated indexing of peaks of the individual domains
- Simultaneous integration of intensities from up to eight individuals, full graphics control
- Intensity scaling based on sets of symmetry-related reflections



## DIAGNOSTIC DIAGRAMS

After scaling intensity data, the improvement can be checked visually



## MULTI-DOMAIN INTEGRATION

Example of a two-domain integration (turquoise: groups of overlapping reflections)

### INCOMMENSURATELY MODULATED STRUCTURES

- Evaluation of the main lattice
- Determination and refinement of up to 3 q vectors
- Integration of main reflections as well as satellites
- Ability to process patterns from multi-domain crystals being incommensurately modulated



SINCE 1887

STOE & Cie GmbH  
Hilpertstrasse 10  
64295 Darmstadt | Germany

P +49 (0) 6151 988 70  
F +49 (0) 6151 988 788  
info@stoe.com

WWW.STOE.COM